

Tissumend II Products

To fit all your tissue adhesive needs

Absorbable methoxypropyl cyanoacrylate composition designed for closing and sealing external tissues:

- Quick 10 second set up
- Strong bond between tissue surfaces
- Hemostatic properties
- Tissue friendly
- Absorption via hydrolysis in 60-90 days
- Can be used with and without sutures
- Easy to use

For external use:

Tissumend II Synthetic Absorbable Tissue Adhesive

Tissumend II in 0.2 ml prefilled pipettes, 20 pipettes per box. Store under refrigeration



tail dock

little lacerations

neuter/spay

dewclaw

paw pad cuts

nail bleeders

hematoma

nose cuts

debrided traumatic wounds

Tissumend II SC Synthetic Absorbable Tissue Adhesive

Tissumend II SC in 1.5 ml multi-use glass bottles with 10 applicator pipettes. No refrigeration required



Tissumend II being applied to close canine castration site




VETERINARY PRODUCTS
LABORATORIES

Innovative Products in Veterinary Medicine™
P.O. Box 34820 • Phoenix, Arizona 85067-4820
www.vpl.com 888-241-9545

Tissumend II Sterile Synthetic Absorbable Tissue Adhesive

The only absorbable tissue adhesive designed to close and seal internal tissues.

For internal use:



Tissumend II Sterile being applied to close lung laceration

Indications:

- Lung
- Liver
- Spleen
- Cornea
- Omentum
- Bronchial tubes
- Pancreas
- Kidney
- Other internal and external surgeries



Deep eroding ulcer with neovascularization in a 3 year old male Chihuahua



After debridement to remove the nonadherent edges, Tissumend II Sterile was applied to the ulcer



Four weeks after application of Tissumend II Sterile

Corneal ulcer photos courtesy of Dr. Rocco Mele, Tucson, AZ



Tissumend II Sterile in 1.2 ml prefilled pipettes, 6 pipettes per box. Store under refrigeration

How Do Surgical Adhesives Hold Together?*

	Tissumend II	Tissumend II SC	Tissumend II Sterile	Nexaband®	Vetbond™
External Use	YES	YES	YES	YES	YES
Internal Use	NO	NO	YES	NO	NO
Sterile	NO	NO	YES	NO	NO
Absorbable	YES	YES	YES	NO	NO
Prefilled Pipettes	YES	NO	YES	NO	NO
Multi-Use Bottles	NO	YES	NO	YES	YES

*Data obtained from Nexaband and Vetbond websites September 2008. Nexaband is registered trademark of Abbott Laboratories. Vetbond is trademark of 3M.

