


PROVEN SEALING FOR FASTER HEALING - INTERNALLY AND EXTERNALLY™

TISSUMEND II STERILE

*The ONLY absorbable tissue adhesive designed
to close and seal internal tissues*



Indications:

- Lung
- Liver
- Spleen
- Cornea
- Bronchial tubes
- Pancreas
- Kidney
- Other internal and external surgeries

Tissumend II Sterile being applied to close lung laceration

Tissumend II Sterile Absorbable Tissue Adhesive is a sterile absorbable methoxypropyl cyanoacrylate composition

- Forms a tough durable bond for internal and external tissues
- Is non-reactive to intact and lacerated tissues
- Enhances healing with no leakage
- Can be used with and without sutures
- Provides quick 10 second set up for most tissues
- Has hemostatic properties
- Absorbs via hydrolysis in 60-90 days
- Is available in easy to use prefilled pipettes to save surgery time



*Can be used with
MicroDose™ syringe*

For more information please call



VETERINARY PRODUCTS
LABORATORIES

Innovative Products in Veterinary Medicine™
P.O. Box 34820 • Phoenix, Arizona 85067-4820

www.vpl.com 888-241-9545

*Now there is one tissue adhesive available for all your surgeries.
NO OTHER TISSUE ADHESIVE can be used both internally and externally.*

TISSUMEND II STERILE

Proven in Studies for Specialized Surgeries*

- The effectiveness of Tissumend II Sterile has been demonstrated in repair of renal lacerations, corneal defects and lung lacerations
- Tissumend II Sterile is a viable alternative for sealing soft tissue defects or injuries instead of suturing

SUMMARY OF STUDIES

Repair of Renal Lacerations

- Tissumend II Sterile forms rapid seal on renal incisions and provides hemostasis
- Tissumend II Sterile helps damaged kidney tissue heal in a timely manner
- Absence of post-operative adhesions makes Tissumend II Sterile a desirable choice for internal soft tissue repair

Repair of Corneal Defects

- Tissumend II Sterile allows approximation of defective cornea tissue and acts as a barrier to allow damaged corneal tissue to heal
- Use of Tissumend II Sterile appears to aid in the restoration of damaged corneal tissue without adverse histological effects

Closure Methods for Lung Lacerations

- Tissumend II Sterile effectively seals lungs without tearing and without leakage, even with active recoveries
- Tissumend II Sterile can be applied over suture lines of tissues under tension, such as tracheal anastomoses, to help prevent postoperative leakage and subcutaneous emphysema

- Specifications**
- Supplied in 6, 0.2 ml pipettes per box
 - Packaged under dry nitrogen for increased stability and shelf life
 - Store under refrigeration



How Do Surgical Adhesives Hold Together?*

	Tissumend II Sterile	Tissumend II	Nexaband®	Vetbond™
Absorbable	YES	YES	NO	NO
Sterile	YES	NO	NO	NO
Internal Use	YES	NO	NO	NO
External Use	YES	YES	YES	YES
Prefilled Pipettes	YES	YES	NO	NO
Syringe Applicator	YES	YES	NO	NO

*Data obtained from Nexaband and Vetbond websites in 03/05. Nexaband® is registered trademark of Abbott Laboratories. Vetbond™ is trademark of 3M.

PROVEN SEALING FOR FASTER HEALING - INTERNALLY AND EXTERNALLY™